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- AKERSTRÖM, B. & BJÖRCK, L. 1986. A physicochemical study of protein G, a molecule with unique immunoglobulin G-binding properties. *Journal of Biological Chemistry*, 261, 10240-10247.
- AKERSTROM, B., NIELSEN, E. & BJORCK, L. 1987. Definition of IgG- and albumin-binding regions of streptococcal protein G. *The Journal of biological chemistry*, 262, 13388-91.
- BJORCK, L. 1988. Protein L. A novel bacterial cell wall protein with affinity for Ig L chains. *Journal of immunology*, 140, 1194-7.
- BJORCK, L. & KRONVALL, G. 1984. Purification and some properties of streptococcal protein G, a novel IgG-binding reagent. *Journal of immunology*, 133, 969-74.
- DEISENHOFER, J. 1981. Crystallographic refinement and atomic models of a human Fc fragment and its complex with fragment B of protein A from *Staphylococcus aureus* at 2.9- and 2.8-Å resolution. *Biochemistry*, 20, 2361-70.
- DERRICK, J. P. & WIGLEY, D. B. 1992. Crystal structure of a streptococcal protein G domain bound to an Fab fragment. *Nature*, 359, 752-4.
- ELIASSON, M., ANDERSSON, R., OLSSON, A., WIGZELL, H. & UHLEN, M. 1989. Differential IgG-binding characteristics of staphylococcal protein A, streptococcal protein G, and a chimeric protein AG. *Journal of immunology*, 142, 575-81.
- ELIASSON, M., OLSSON, A., PALMCRANTZ, E., WIBERG, K., INGANAS, M., GUSS, B., LINDBERG, M. & UHLEN, M. 1988. Chimeric IgG-binding receptors engineered from staphylococcal protein A and streptococcal protein G. *The Journal of biological chemistry*, 263, 4323-7.
- GOMEZ, M. I., LEE, A., REDDY, B., MUIR, A., SOONG, G., PITT, A., CHEUNG, A. & PRINCE, A. 2004. *Staphylococcus aureus* protein A induces airway epithelial inflammatory responses by activating TNFR1. *Nature medicine*, 10, 842-8.
- GOMEZ, M. I., O'SEAGHDHA, M., MAGARGEE, M., FOSTER, T. J. & PRINCE, A. S. 2006. *Staphylococcus aureus* protein A activates TNFR1 signaling through conserved IgG binding domains. *The Journal of biological chemistry*, 281, 20190-6.
- GOUDA, H., TORIGOE, H., SAITO, A., SATO, M., ARATA, Y. & SHIMADA, I. 1992. Three-dimensional solution structure of the B domain of staphylococcal protein A: comparisons of the solution and crystal structures. *Biochemistry*, 31, 9665-9672.
- GOWARD, C. R., MURPHY, J. P., ATKINSON, T. & BARSTOW, D. A. 1990. Expression and purification of a truncated recombinant streptococcal protein G. *The Biochemical journal*, 267, 171-7.
- GRAILLE, M., STURA, E. A., CORPER, A. L., SUTTON, B. J., TAUSSIG, M. J., CHARBONNIER, J. B. & SILVERMAN, G. J. 2000. Crystal structure of a *Staphylococcus aureus* protein A domain complexed with the Fab fragment of a human IgM antibody: structural basis for recognition of B-cell receptors and superantigen activity. *Proceedings of the National Academy of Sciences of the United States of America*, 97, 5399-404.
- GRONENBORN, A. M., FILPULA, D. R., ESSIG, N. Z., ACHARI, A., WHITLOW, M., WINGFIELD, P. T. & CLORE, G. M. 1991. A novel, highly stable fold of the immunoglobulin binding domain of streptococcal protein G. *Science*, 253, 657-61.
- GULICH, S., UHLEN, M. & HOBER, S. 2000. Protein engineering of an IgG-binding domain allows milder elution conditions during affinity chromatography. *Journal of biotechnology*, 76, 233-44.
- GUSS, B., ELIASSON, M., OLSSON, A., UHLEN, M., FREJ, A. K., JORNVALL, H., FLOCK, J. I. & LINDBERG, M. 1986. Structure of the IgG-binding regions of streptococcal protein G. *The EMBO journal*, 5, 1567-75.
- HARTLEIB, J., KOHLER, N., DICKINSON, R. B., CHHATWAL, G. S., SIXMA, J. J., HARTFORD, O. M., FOSTER, T. J., PETERS, G., KEHREL, B. E. & HERRMANN, M. 2000. Protein A is the von Willebrand factor binding protein on *Staphylococcus aureus*. *Blood*, 96, 2149-56.
- HUSE, K., BOHME, H. J. & SCHOLZ, G. H. 2002. Purification of antibodies by affinity chromatography. *Journal of biochemical and biophysical methods*, 51, 217-31.
- INGANAS, M. 1981. Comparison of mechanisms of interaction between protein A from *Staphylococcus aureus* and human monoclonal IgG, IgA and IgM in relation to the classical FC gamma and the alternative F(ab')₂ epsilon protein A interactions. *Scandinavian journal of immunology*, 13, 343-52.
- INGANAS, M., JOHANSSON, S. G. & BENNICH, H. H. 1980. Interaction of human polyclonal IgE and IgG from different species with protein A from *Staphylococcus aureus*: demonstration of protein-A-reactive sites located in the Fab'2 fragment of human IgG. *Scandinavian journal of immunology*, 12, 23-31.
- JANSSON, B., UHLÉN, M. & NYGREN, P.-Å. 1998. All individual domains of staphylococcal protein A show Fab binding. *FEMS Immunology & Medical Microbiology*, 20, 69-78.
- JONES, J. V. 1990. Staphylococcal Protein A as an extracorporeal immunosorbent: theoretical and practical considerations. *Transfusion science*, 11, 153-159.
- KASTERN, W., HOLST, E., NIELSEN, E., SJÖBRING, U. & BJÖRCK, L. 1990. Protein L, a bacterial immunoglobulin-binding protein and possible virulence determinant. *Infection and immunity*, 58, 1217-1222.
- KASTERN, W., SJOBRING, U. & BJORCK, L. 1992. Structure of peptostreptococcal protein L and identification of a repeated immunoglobulin light chain-binding domain. *The Journal of biological chemistry*, 267, 12820-5.

- KRAULIS, P. J., JONASSON, P., NYGREN, P.-Å., UHLÉN, M., JENDEBERG, L., NILSSON, B. & KÖRDEL, J. 1996. The serum albumin-binding domain of streptococcal protein G is a three-helical bundle: a heteronuclear NMR study. *FEBS letters*, 378, 190-194.
- LANGONE, J. J. 1982. Applications of immobilized protein A in immunochemical techniques. *Journal of immunological methods*, 55, 277-96.
- MARKELA, E., STAHLBERG, T. H. & HEMMILA, I. 1993. Europium-labelled recombinant protein G. A fast and sensitive universal immunoreagent for time-resolved immunofluorometry. *Journal of immunological methods*, 161, 1-6.
- MERINO, N., TOLEDO-ARANA, A., VERGARA-IRIGARAY, M., VALLE, J., SOLANO, C., CALVO, E., LOPEZ, J. A., FOSTER, T. J., PENADES, J. R. & LASA, I. 2009. Protein A-mediated multicellular behavior in *Staphylococcus aureus*. *Journal of bacteriology*, 191, 832-43.
- MONGODIN, E., BAJOLET, O., HINNRASKY, J., PUCHELLE, E. & DE BENTZMANN, S. 2000. Cell wall-associated protein A as a tool for immunolocalization of *Staphylococcus aureus* in infected human airway epithelium. *The journal of histochemistry and cytochemistry : official journal of the Histochemistry Society*, 48, 523-34.
- MURPHY, J. P., DUGGLEBY, C. J., ATKINSON, M. A., TROWERN, A. R., ATKINSON, T. & GOWARD, C. R. 1994. The functional units of a peptostreptococcal protein L. *Molecular microbiology*, 12, 911-20.
- NILSON, B., BJORCK, L. & AKERSTROM, B. 1988. Enzyme linked immunosorbent assay using alkaline phosphatase conjugated with streptococcal protein G. *Journal of immunoassay*, 9, 207-25.
- NILSON, B. H., LOGDBERG, L., KASTERN, W., BJORCK, L. & AKERSTROM, B. 1993. Purification of antibodies using protein L-binding framework structures in the light chain variable domain. *Journal of immunological methods*, 164, 33-40.
- NILSSON, B., MOKS, T., JANSSON, B., ABRAHMSSEN, L., ELMLAD, A., HOLMGREN, E., HENRICHSON, C., JONES, T. A. & UHLEN, M. 1987. A synthetic IgG-binding domain based on staphylococcal protein A. *Protein engineering*, 1, 107-13.
- O'SEAGHDHA, M., VAN SCHOOTEN, C. J., KERRIGAN, S. W., EMSLEY, J., SILVERMAN, G. J., COX, D., LENTING, P. J. & FOSTER, T. J. 2006. *Staphylococcus aureus* protein A binding to von Willebrand factor A1 domain is mediated by conserved IgG binding regions. *The FEBS journal*, 273, 4831-41.
- OLSSON, A., ELIASSON, M., GUSS, B., NILSSON, B., HELLMAN, U., LINDBERG, M. & UHLEN, M. 1987. Structure and evolution of the repetitive gene encoding streptococcal protein G. *European journal of biochemistry / FEBS*, 168, 319-24.
- PATELLA, V., CASOLARO, V., BJORCK, L. & MARONE, G. 1990. Protein L. A bacterial Ig-binding protein that activates human basophils and mast cells. *Journal of immunology*, 145, 3054-61.
- SAUER-ERIKSSON, A. E., KLEYWEGT, G. J., UHLEN, M. & JONES, T. A. 1995. Crystal structure of the C2 fragment of streptococcal protein G in complex with the Fc domain of human IgG. *Structure*, 3, 265-78.
- SJOBRING, U., BJORCK, L. & KASTERN, W. 1991. Streptococcal protein G. Gene structure and protein binding properties. *The Journal of biological chemistry*, 266, 399-405.
- SJODAHL, J. 1977. Structural studies on the four repetitive Fc-binding regions in protein A from *Staphylococcus aureus*. *European journal of biochemistry / FEBS*, 78, 471-90.
- UHLEN, M., GUSS, B., NILSSON, B., GATENBECK, S., PHILIPSON, L. & LINDBERG, M. 1984. Complete sequence of the staphylococcal gene encoding protein A. A gene evolved through multiple duplications. *The Journal of biological chemistry*, 259, 1695-702.
- WIKSTROM, M., SJOBRING, U., KASTERN, W., BJORCK, L., DRAKENBERG, T. & FORSEN, S. 1993. Proton nuclear magnetic resonance sequential assignments and secondary structure of an immunoglobulin light chain-binding domain of protein L. *Biochemistry*, 32, 3381-6.
- YOUNG, W. W., JR., TAMURA, Y., WOLOCK, D. M. & FOX, J. W. 1984. Staphylococcal protein A binding to the Fab fragments of mouse monoclonal antibodies. *Journal of immunology*, 133, 3163-6.